

Compact
Powerful
Precise



ColorFlex[®] L2 Coffee

*Level UP your Color Confidence with the new
ColorFlex L2 Coffee Spectrophotometer.*

As simple as you want, as powerful as you need.

Elevate your color accuracy and user experience like never before. Compact and ergonomic build requires minimal space while offering enhanced durability and reliability. Designed to streamline and expedite your quality control processes, ColorFlex L2 Coffee Spectrophotometer is expertly crafted to measure the color of ground coffee.




Explore the Next-Gen ColorFlex®...L2 Coffee



#1 Color Measurement Solution for Gound Coffee

HunterLab Coffee Reference Standard: Coffee tile to be used as a Performance Qualification step with values Read-at Factory in HCCI and SCCA- Hunter Coffee Color Index; SCAA #, SCAA Roast Classification and CIE L*, a*, b* values.




All the Coffee Indices You Need

Measure Coffee in Every Form: HCCI (HunterLab Coffee Color Index values), SCAA Number, SCAA Roast Classification.

45°a:0° Annular Illuminated Optical Geometry

Superior Consistency and Color Accuracy: Enhancing the efficiency of Quality Control workflows and delivery of superior quality products repeatedly.



Compact with Sealed Spillproof Case

Space-Saving and Easy Maintenance: Occupies minimal space and is designed to withstand the demands of busy lab environments.



COFFEE



Powerful Standalone Functionality

No computer required: With huge onboard processing and storage capacity, ColorFlex L2 Coffee can function as an independent color workstation. Connect monitor, keyboard, and mouse if desired.




Modern Tablet User Interface with Touchscreen Color Display

Durable and Efficient: Facilitates accurate color measurement in seconds, with a robust industrial glass design suitable for harsh environments.




New & Improved Color Software

Essentials 2.0: Take precise and reliable color measurement in seconds with next-gen user-friendliness and seamlessly analyze multiple tristimulus and spectral data sets. It's that powerful.




Self-Guided and Intuitive Out-of-Box Startup

Time-Saving and User-Friendly: Ready to use straight out of the box with an integrated wizard for quick training, reducing operator time and easing concerns related to staff training and turnover.



Built-in Advanced Diagnostics

On-Board Diagnostic Health Check: Ensures the HunterLab ColorFlex L2 Coffee is always operating at peak performance, continuing to build on the legendary history of reliability of the ColorFlex family of instruments.



COFFEE



Advanced Data Export and Communication

Full suite of data communication options: Features USB ports, HDMI, and support for mouse and keyboard, enabling easy data export/printing.

RJ45 Connectivity for LIMS, SPC, and Data Collection

Simplifies QC workflows: Aids in trend analysis and workflow simplification across manufacturing sites.



CONTACT US

We've worked with thousands of companies, but you'll feel like our one and only. If color accuracy is crucial to your business, you've come to the right place—with the right people. Let's get in touch. Ask a question, request a quote or just say hello!

OUR OFFICES & DISTRIBUTION CENTERS



HunterLab

11491 Sunset Hills Road
Reston, VA 20190-5280
Website: www.hunterlab.com
Email: sales@hunterlab.com
Ph: 703-471-6870
Fax: 703-471-4237

Corporate Headquarters



HunterLab Europe

Christian Jansen, Managing Director
Dr.-August-Einsele-Ring 15
82418 Murnau, Germany
Website: www.hunterlab.com
Email: Christian.Jansen@hunterlab.com
Ph: +49 8841 9464

European Office



HunterLab Asia-Pacific

Heng Su Kiang, APAC Sales Manager
28-2, Plaza KLTS
No. 99, Jalan Gombak
Setapak
53000 Wilayah Persekutuan Kuala Lumpur
Malaysia
Website: www.hunterlab.com
Email: sales-asiapacific@hunterlab.com
Ph: +(60) 340 324 749



HunterLab Asia-Pacific

Leo Cheung, Regional Manager
Flat E, 13F Leahander Centre
28 Wang Wo Tsai Street Tsuen Wan
New Territories Hong Kong
Website: www.hunterlab.com
Email: sales-asiapacific@hunterlab.com
Ph: +852-2406 1723

Asia-Pacific Offices

For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com

HunterLab and ColorFlex are Trademarks of Hunter Associates Laboratory, Inc.
Specifications subject to change without notice

ISO 17025 Certified; ISO 9001 Certified; Certified